

## KIT FOR THE DETECTION OF M235T POLYMORPHISM OF THE ANGIOTENSINOGEN (AGT) GENE

**AMPLI-SET-AGT M235T**

**Cat. n. 1.354**

Products encoded by genes involved in renine-angiotensin-aldosterone system (RAAS) are the natural candidates for the maintenance of sodium homeostasis and regulation of blood pressure. The RAAS system has been shown to be involved in many cardiovascular diseases, including myocardial fibrosis and hypertrophy in hypertensive heart disease, congestive heart failure, myocardial infarction and cardiomyopathy.

Many polymorphisms have been detected in genes of the renine-angiotensin-aldosterone system, as insertion/deletion (I/D) of the ACE gene, the polymorphisms G-217A, G-152A, A-20C, G-6A, M235T and T174M of the angiotensinogen (AGT) gene and the A1166C polymorphism of the angiotensin II type I receptor gene.

The kit allows the detection of the M235T of the AGT gene. The detection of the polymorphism is performed with the amplification with specific primers of a fragment of 303 bp, followed by restriction section due to Lwe I. The polymorphism of a single nucleotide (SNP) in exon 2 causes the substitution of methionine with threonine. In a control population, the frequencies of allele M is 64%, whereas the frequencies of allele T is 36%. In many studies the variant TT has been associated to development of hypertension. Moreover the study of the association of polymorphisms of many genes involved in the RAAS system is very important.

**Principle of Assay:** A) extraction of genomic DNA B) amplification C) enzymatic digestion D) detection on agarose gel.

**Applicability:** On extracted and purified genomic DNA from whole blood samples.

**Tests:** 45.

### ANALYSIS OF RESULTS

The amplified product is a fragment of 303 bp; digested product of the normal allele MM gives two fragments of 266 bp and 37 bp, whereas the mutant allele TT gives a fragment of 303 bp.

### REAGENTS AND STORAGE

#### AMPLIFICATION and DIGESTION

PCR mix	-20°C
H <sub>2</sub> O sterile	-20°C
Taq Polymerase (5U/μl)	-20°C
Lwe I enzyme (10U/μl)	-20°C
Digestion buffer 10X	-20°C
Positive Heterozygous control	-20°C

**Stability:** over 12 months if correctly stored.

1 Wild Type subject MM	2 Heterozygous subject MT	3 Homozygous mutant subject TT
2 bands	3 bands	1 band
266 bp 37 bp	303 bp 266 bp 37 bp	303 bp

#### References:

*Circ Res.* 1990; 67:1355-1364.  
*J Clin Invest.* 1993; 91: 1269-1281.  
*Lancet* 1994; 343:125.  
*J Hypertens.* 1995; 13:1602-1609.  
*J Clin Invest.* 1997;99:1786-1797.  
*Circulation* 2004; 109: 1640-1646.  
*Experimental and Molecular Medicine* 2003; 35, 6, 545-549.  
*Hypertension* 2004; 43:598-602.  
*BMC Cardiovascular Disorders* 2005, 5:7.  
*Obesity Research* 2004, 12, 3, 556-561.